
Sugi®

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

Sugi®
nasal surgery swabs with retrieval cord
according to Prof. Dr. H. Stammberger

12 x 10 mm

100 swabs (10 x 10 pieces, in sachets) REF 31731	60 swabs (10 x 6 pieces, in sachets) REF 31732
---	---



15 x 13 mm

100 swabs (10 x 10 pieces, in sachets) REF 31741	60 swabs (10 x 6 pieces, in sachets) REF 31742
---	---

Sugi® Rhino-Long
double layers with retrieval cord (200 mm)
40 x 12 mm
 40 swabs
 (10 x 4 pieces, in sachets)
 REF 31778

Sugi® Endopads
Ø 34 mm
according to Prof. Dr. H. Stammberger
 (6 x 1 piece, in sachets)
 REF 31785

Sugi® OtoSwab
with retrieval cord
Ø 3 x 1 mm, expanded 3 x 5 mm
 60 swabs (10 x 6 pieces, in sachets)
 REF 31651

Composition

Cellulose: approx. 35% (m/m)
 Regenerated cellulose: approx. 65% (m/m)
 Cord: polyester

Properties

- Sugi absorbent materials are manufactured from high-quality cellulose material and are toxicologically harmless.
- The swabs/pads expand immediately on contact with liquid and become soft and pliable.
- Sugi absorbent materials exhibit high (> 50 N) tear strength.
- Sugi absorbent materials do not adhere to the wound surface and are therefore easily removed.
- Sugi absorbent materials leave a clean wound surface after removal.
- The absorptive capacity is approximately 17 times its own weight.

Range of Applications

a) Nasal surgery swabs

- Preoperative and intraoperative nasal swabs for use in nasal surgery especially during:
 - Septoplasty
 - Septorhinoplasty
 - Conchotomy
 - Paranasal sinus surgery
- Swabs for use in endoscopic ENT surgery
- Swabs for use in skull base surgery involving nasal surgery procedures
- Carrier medium for liquid antibiotic or vasoconstrictor medication

b) Endopads

- Intraoperative cleaning of optical instruments (e.g. endoscopes) and simultaneous application of an anti-fogging agent.

c) OtoSwab

- Intra-operative swab for ear surgery
 - Post-auricular approach to middle ear surgery
 - Endaural approach to middle ear surgery
 - Meatoplasty
 - Tympanoplasty
 - Tympanotomy
 - Cortical Mastoidectomy
 - Middle ear surgical procedures
 - Cochlear Implant
- Carrier medium for vasoconstrictive/antibiotic solutions

Contraindications and Incompatibility

None known

Application

a) Nasal surgery swabs */OtoSwab

The swabs are soaked in a vasoconstrictor and/or anaesthetic solution for preparing the operation site. They are then placed on the relevant mucosal areas for 5 – 10 minutes.

They are easily removed using the attached retrieval cord.

The swabs can also be used intraoperatively to absorb blood and secretion.

b) Endopads

The dry pad is prepared by soaking it in an anti-fogging agent. After removing the protective foil from the adhesive strip, the pad is positioned close to the operation site.

Note

- Sugi absorbent materials should be disposed of like bandaging materials.
- Sterility is only guaranteed if the packaging is undamaged.
- Storage: store dry at room temperature (18 - 25 °C).
- The shown Sugi absorbent materials are radiolucent.
- Sugi absorbent materials are for use by qualified medical personnel only.
- * In rare instances the physiochemical conditions associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include: sudden fever (usually 102° or more), vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn. Should one of these symptoms occur, administer the appropriate medical treatment.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

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